

# Quadrant<sup>®</sup> LSG PPSU – Life Science Grade (LSG) Material Solutions for Medical Applications

# Plastic stock shapes biocompatible according to USP and ISO 10993-1

### Trends

In addition to extreme resistance to steam and chemical sterilisation, today's medical device designs require exceptional physical performance that can help deliver the extended product life to meet conscious users demand.

## **Quadrant Answers**

Quadrant<sup>®</sup> LSG PPSU is an outstandingly tough engineering plastic with a temperature of deflection under load of about 205 ° C, excellent chemical resistance, dimensional stability and the ability to withstand more than 1,000 cycles of steam sterilisation without significant loss of properties. It is highly resistant to cleaning agents, disinfectants and various solvents, as well as to gamma radiation.

#### **Customer Benefits**

Quadrant delivers unique LSG PPSU materials for machining which were tested on stock shapes in order to check their compliance with requirements for biocompatibility of materials per United States Pharmacopeia (USP) and ISO 10993-1. This offers the designer increased safety, approval facilitation, as well as reduced cost and time to market. Many standard colours including blue, grey,natural, black, red, green, orange, rust and yellow are offered as rod and plate. Various combinations of shape sizes provide manufacturers an ideal yield with minimal waste.

## Properties of Quadrant® LSG PPSU by Quadrant

- Superior resistance to common sterilisation techniques compared to other engineering plastics
- Extremely high resistance to repeated steam sterilisation
- High resistance to gamma radiation
- Tested biocompatibility per USP and ISO 10993-1 on stock shapes
- Excellent dimensional stability
- Very high impact strength
- Excellent machinability
- Lot traceability



### **Typical applications**

- Trial implants
- Medical instruments
- Endoscopic equipment
- Radiation therapy and brachytherapy components
- Trays for trauma products
- Biotechnology and laboratory equipment
- Analytical and diagnostic equipment

#### **Biocompatibility status of Quadrant® LSG PPSU grades:**

TESTS (1) (2) MATERIALS	1. Cytotoxicity Bef. ISO 10993-5 and USP <87> Biological	tivity, In Vitro Elution Test ensitization		3. Intracutaneous Reactivity	Ref.: ISO 10993-10 and USP <88>, Biological Reactivity Tests, In Vivo - Intracutaneous Test	4. Acute Systemic Toxicity	Ref.: ISO 10993-11 and USP <88>, Biological Reactivity Tests, In Vivo - Systemic Injection Test	5. Implantation Test	Ref.: USP <88>, Biological Reactivity Tests, In Vivo - Implantation Test (7 days)	6. Human blood compability	Ref.: ISO 10993-4, Indirect Hemolysis (in vitro)	7. USP-Physicochemical Tests for Plastic	Ref.: USP <661 > Containers, Ultra Pure Water extract 70°C/24h	USP Class VI (conclusion from tests 3,4 and 5)	
Quadrant <sup>®</sup> LSG PPSU black	$\checkmark$		$\checkmark$		$\checkmark$		$\checkmark$		$\checkmark$		$\checkmark$		$\checkmark$	$\checkmark$	
Quadrant <sup>®</sup> LSG PPSU natural (ivory)	$\checkmark$		V(3)		$\checkmark$		$\checkmark$		V(3)		V(3)		$\checkmark$	V(3)	
Quadrant <sup>®</sup> LSG PPSU blue, green, grey, orange, red, rust, yellow	$\checkmark$		V(3)		V(3)		V(3)		V(3)		V(3)		$\checkmark$	V(3)	

✓ This test was carried out and the material passed the test. V(3) An abbreviated biocompatibility type testing programme was run on these Quadrant® LSG PPSU products, which are colour variations of Quadrant® LSG PPSU black. A comprehensive biocompatibility type testing programme was run on Quadrant® LSG PPSU black. (2) Quadrant EPP performs testing on its Life Science Grades in order to facilitate evaluation by its customers of their biocompatibility with regard to the requirements applicable to the specific use of the finished product. Quadrant EPP does not possess expertise in evaluating the suitability of its tested materials for use in specific medical, pharmaceutical, or biotechnological applications. It remains the customer's sole responsibility to test and assess the suitability of Quadrant's Life Science Grades for its intended applications, processes and uses.

#### All tests were run on test specimens machined from rod diameter 50 mm shortly after manufacture.

## **Quadrant Engineering Plastic Products**

#### Europe

Quadrant EPP Deutschland GmbH Max-Planck-Straße 11 48691 Vreden, Germany T +49[0] 2564 3010 F +49[0] 2564 3012 55 contact@qplas.com North America Quadrant EPP USA, Inc. 2120 Fairmont Avenue PO Box 14235 - Reading, PA 19612-4235 T 800 366 0300 | +1 610 320 6600 F 800 366 0301 | +1 610 320 6638 americas.epp@qplas.com

#### Asia-Pacific

Available in colours: natural (ivory), black, blue, green, grey, orange,

red, rust, yellow, as rod

and plate, and also as

finished parts. Please

consult us for availa-

bility of rod diameters

and plate sizes.

Quadrant EPP Asia Pacific Ltd 60 Ha Mei San Tsuen, Ping Shan Yuen Long - NT Hong Kong T +852 24702683 F +852 24789966 asia.epp@qplas.com

This brochure and any data and specifications presented here or on our website shall provide promotional and general information about the Engineering Plastic Products (the "Products") manufactured and offered by Quadrant Engineering Plastic Products ("Quadrant") and shall serve as a preliminary guide. All data and descriptions relating to the Products are of a general informational nature only. Neither this brochure nor any data and specifications presented on our website shall create or be implied to create any legal or contractual obligation. This brochure and any dato are specifications herein do not create expressly or by implication any legal, contractual or warranty obligation whatsoever. No warranty of any kind, either express or implied, is made as to the information contained in these pages, including, but not limited to, all warranties provided for by Louisian law, any implied warranty of merchantability, of fitness for a particular purpose, and any warranty against hidden defects or redhibitory defects or vices. No information in this brochure creates any express or implied warranty that the goods described herein shall conform to any description herein. Quadrant sells the products described herein syle to sophisticated users and not to consumers, and Quadrant assumes no responsibility that any goods described herein unter contract.

Any illustration of the possible fields of application of the Products shall merely demonstrate the potential of these Products, but any such description does not constitute any kind of covenant or warranty whatsoever. Irrespective of any tests that Quadrant may have carried out with respect to any Product, Quadrant does not possess expertise in evaluating the suitability of its materials or Products for use in specific applications or products manufactured or offered by the customer respectively. It thus remains the customer's sole responsibility to test and assess the suitability and compatibility of Quadrant's Products for its intended applications, processes and uses, and to choose those Products that according to its assessment meet the requirements applicable to the specific use of the finished product. The customer undertakes all liability in respect of the application, processing or use of the aforementioned information or product, or any consequence thereof, and shall verify its quality and other properties.

Quadrant's Products should not be used for applications involving medical devices that are intended to remain implanted in the human body continuously for a period exceeding 24 hours (30 days"), or are intended to remain in contact with internal human tissue or bodily fluids for more than 24 hours (30 days"), or as critical components of medical devices that are essential to the continuation of human life.

components of medical devices that are essential to the continuation of numan line. Quadrant is not a medical device manufacturer and the information herein does not constitute any express or implied warranties or representations whatsoever, including, but not limited to, all warranties provided for by Louisiana law, any implied warranty of merchantability, of fitness for a particular purpose, any warranty against hidden defects or redhibitory defects or vices, or that Quadrant's materials are manufacturered in accordance with the quality standards appropriate and necessary for materials intended for use in implantable medical device applications and in applications that are essential to the restoration of or continuation of a bodily function important to the continuation of human life. \*: "30 days" applies to Ketron® PEEK-CLASSIXTM LSG white only.

Quadrant® is a registered trademark of the Quadrant Group of companies.

# More information: www.quadrantplastics.com



Lit.LSG\_PPSU\_04/2014

# You inspire ... we materialize®